




KN95 FACE MASK



Por favor, hagan sus pedidos enviando un correo electrónico a msamper@compradorasesor.es


Main Material:

1. Outer 2 layer non-woven fabric, filter large particle dust and fuzz.
2. Middle 2 layer melt blown non-woven fabric, more than 95% particles filtration efficiency.

Strong filtration, adsorb PM2.5, haze, bacteria and virus..

3. Inner: hot air fabric, bacteriostatic material, eliminate secondary pollution.

Product PFE \geq 95%
Reached KN95 level



10pcs/polybag, 20 polybags / color box, 50box/carton。

Total: 1000pcs/carton

G.W:8.4Kg, N.W:7.4Kg

Carton size: 62cm*33cm*47cm, 50box/carton



FACE MASK

EFFECTIVE PROTECTION OF MICRO PARTICLE DUST

Usage Method:

1. Keep the skin dry, open the mask, the light-colored cloth toward the face and the nose bar above;
2. Hang the ear sling on the left and right sides of both ears and adjust, so that the force on both ears is even;
3. Adjust the position of the mask, spread the mask up and down, completely cover the mouth and nose.

Product Description:

- High quality polypropylene non-woven fabric, 3 layer thickening protection design
- Low sensitive, non-irritation material
- Low resistance to breathing
- Nose bar adaptable
- High filtration capacity
- Perfect fitting

Wenzhou

Co.,Ltd

Material: Non woven fabric, Melt blown fabric, Not an cotton

Standard: GB 2626-2006

Date of Manufacture:

Date of Expiry: 2 years

Attention:

- Please wear the mask according to the picture instructions and check the tightness between the mask and the face;
- If the mask is damaged, please do not continue to use it. It is recommended to replace it immediately.
- Do not use the product at bedtime.
- After wearing the mask, please avoid approaching the fire source or heat source, which may cause burning or deformation of the material.
- If it needs to be removed during use, please avoid the external contact of the mask with the nose and mouth.
- This product is disposable mask, do not reuse it.
- After use, please fold the mask from inside to outside and tie it with ear strap before putting it into the garbage bin.
- It is not recommended to be used in places with strict microbial index control.



Por favor, hagan sus pedidos enviando un correo electrónico a msamper@compradorasesor.es


专业生产

KN95口罩车间生产场景



专业生产

KN95口罩车间生产场景



Por favor, hagan sus pedidos enviando un correo electrónico a msamper@compradorasesor.es



Por favor, hagan sus pedidos enviando un correo electrónico a msamper@compradorasesor.es

China mask standard

(GB-32610-2016 / GB2626-2006)

Certification/ Class (Standard)	N95 (NIOSH-42C FR84)	FFP2 (EN 149-2001)	KN95 (GB2626-20 06)	P2 (AS/NZ 1716:2012)	Korea 1 st Class (KMOEL - 2017-64)	DS (Japan JMWHLW- Notification 214, 2018)
Filter performance – (must be ≥ X% efficient)	≥ 95%	≥ 94%	≥ 95%	≥ 94%	≥ 94%	≥ 95%
Test agent	NaCl	NaCl and paraffin oil	NaCl	NaCl	NaCl and paraffin oil	NaCl
Flow rate	85 L/min	95 L/min	85 L/min	95 L/min	95 L/min	85 L/min
Total inward leakage (TIL)* – tested on human subjects each performing exercises	N/A	≤ 8% leakage (arithmetic mean)	≤ 8% leakage (arithmetic mean)	≤ 8% leakage (individual and arithmetic mean)	≤ 8% leakage (arithmetic mean)	Inward Leakage measured and included in User Instructions
Inhalation resistance – max pressure drop	≤ 343 Pa	≤ 70 Pa (at 30 L/min) ≤ 240 Pa (at 95 L/min) ≤ 500 Pa (clogging)	≤ 350 Pa	≤ 70 Pa (at 30 L/min) ≤ 240 Pa (at 95 L/min)	≤ 70 Pa (at 30 L/min) ≤ 240 Pa (at 95 L/min)	≤ 70 Pa (w/valve) ≤ 50 Pa (no valve)
Flow rate	85 L/min	Varied – see above	85 L/min	Varied – see above	Varied – see above	40 L/min
Exhalation resistance - max pressure drop	≤ 245 Pa	≤ 300 Pa	≤ 250 Pa	≤ 120 Pa	≤ 300 Pa	≤ 70 Pa (w/valve) ≤ 50 Pa (no valve)
Flow rate	85 L/min	160 L/min	85 L/min	85 L/min	160 L/min	40 L/min
Exhalation valve leakage requirement	Leak rate ≤ 30 mL/min	N/A	Depressurizatio n to 0 Pa ≥ 20 sec	Leak rate ≤ 30 mL/min	visual inspection after 300 L/min for 30 sec	Depressurizatio n to 0 Pa ≥ 15 sec
Force applied	-245 Pa	N/A	-1180 Pa	-250 Pa	N/A	-1,470 Pa
CO ₂ clearance requirement	N/A	≤ 1%	≤ 1%	≤ 1%	≤ 1%	≤ 1%

*Japan JMWHLW-Notification 214 requires an Inward Leakage test rather than a TIL test.



1 B0011112242



(2018) 国认监认字(244)号



检 验 报 告

TEST REPORT



扫一扫关注我们



扫一扫查询真伪



报 告 编 号

REPORT NO

Q @ @ W202005430 @

产 品 名 称

KN95 @

NAME OF SAMPLE

委 托 单 位

C U S T O M E R

检 验 类 别

委托检验

TEST CATEGORY

浙江省轻工业品质量检验研究院

(浙江省纺织测试研究院)

Zhejiang Light Industrial Products Inspection and Research Institute

National Textiles and Garment Quality Supervision Inspection Center(Zhejiang)

检验报

JJ W2O200S430 Q

第 1 页 共 3 页

Name of Customer		Address	
J4anufac turer		Address	
pp iq Sai0ple informal ion	@iPnQJ Name of sainple: 095 fi @iPn@J Charac ter ist i es: Q & Trademark: --- H /Q H Spec if icat ion/node1: --- Level: 095 f Category of safety speci f i cat i on: 样品款号/货号 Art. No.: ---		
	t _LJ Q ,N (ft bove-v ent i oned informat ion by Customer-supplied)		
The sent hay of sai0ple	自送	样品数量 Sample quant ity	19 @
Recei ving Oate ol Sainple	2020-0410	检测类别 Test Category	委托检验
判定依据 RatingRequirements	GB 2626—2ü06		
@áf@Ö/Test Suninary:			
(检验报告专用章) Test Seal			
}Íuf B J/ Oate of Approval: 2020-04-12			
Remarks			

Approved' by

报

JJ W202005430 Q

第 2 页 共 3 页

		检测方法	单位	标准要求 (KN95)	实测值	单项评价	结果备注
1	@@@ @ @	GB 2626-2006 6.3	8	195.0	95.0	符合	---
2	呼吸阻力	GB 2626-2006 6.6	Pa	T350	93	符合	---
	呼气阻力		Pa	250	66		
3	外观检查	GB 2626-2006 6.1 6.2		符合 GB2626-2006 标准要求	2 个口罩的外 观质量符合标 准	符合	
4	O @T @@	GB 2890-2009 6.8		N 60	s	符合	

检验报

@ @ @ W202005430 @

第 3 页 共 3 页



声明

一、本机构保证检测的公正性、独立性和诚实性，对检测结果负责，对委托方所提供的检测样品保密和保护所有权。
二、本报告无签发人员签字或未加盖本机构红色检验检测专用章无效。
三、报告涂改无效。复印件未重新加盖本机构红色检验检测专用章无效。
四、本报告的检测数据和结果仅对送检样品负责。
五、本报告各页均为报告不可分割之部分，使用者单独抽出某些页导致误解或用于其它用途及由此造成的后果，本机构不负相应的法律责任。
六、委托方若对本报告有异议，应及时向本机构提出。政府行政管理部门下达的指令性任务，被检方对检验结果有异议时，应按政府行政管理部门文件规定及国家相关法律、法规规定进行。

DECLARATION

1. Our organization guarantees impartiality, independence and honesty of inspection, and is responsible for the results of inspection, keeping the samples supplied by the entrusting party confidential and at the same time protecting the ownership of the samples supplied.
2. The test report will be deemed invalid without signatures of the inspector/reviewer and authorized personnel, and the red special inspection stamp of our organization.
3. The test report will be invalid if it is altered. Copies of the report are invalid without the red special inspection stamp of our organization.
4. The test results shown in this report is only valid for the tested sample.
5. All the pages of the report are integral parts of the report. Our organization will not be responsible for any misunderstanding or other results caused by using separate page(s) of the report.
6. If there is any dissent of the report, the entrusting party shall notify our organization timely. For the mandatory inspection given by governmental administration departments, and dissent about the sample being tested or test results on the report should be dealt with in accordance with national regulations.

The Affiliated Inspection Centers

浙江省轻工业产品质量检验研究院
Zhejiang Light Industrial Products Inspection and Research Institute
国家纺织品和服装质量监督检验中心(浙江)
National Textiles and Garment Quality Supervision Inspection Center (Zhejiang)
地址: 浙江省杭州市下沙路300号6楼
Address: Building No.6, 300 XiaSha Road, Hangzhou, Zhejiang

电话: 0571-85122669
Telephone: 0571-85122669
E-Mail: 1137907994@qq.com

国家家具产品质量监督检验中心(浙江)
National Center for Quality Supervision Inspection of Furniture (Zhejiang)
浙江室内安全及家具产品质量检验中心
& 浙江室内安全及家具产品质量检验中心
地址: 浙江省宁波市江北区江东南路4号
Address: Building 4 LiangZhu University Science and Technology Part No.1, JingYi Rd, Yuhang District, Hangzhou, Zhejiang

电话: 0571-89009556
Telephone: 0571-89009556
E-Mail: 2047699564@qq.com

浙江省轻工及五金产品质量检验中心
Zhejiang Center of Quality Test for Light Industry and Hardware Products
浙江体育用品质量监督检验中心
SU : QIT Jili iji UN X & m 222 3
地址: 浙江省宁波市江北区江东南路222号
Address: No. 222 Tianmushan Rd., Hangzhou, Zhejiang
Building 3 LiangZhu University Science and Technology Park, No.1, JingYi Rd., Hangzhou, Zhejiang

电话: 0571-89001107
Telephone: 0571-89001107
E-Mail: vnijtest@126.com

国家锁具产品质量监督检验中心(浙江)
National Center for Quality Supervision Inspection of Lock (Zhejiang)
浙江锁具产品质量检验中心
SU vidi UI- N 2s U
地址: 浙江省宁波市江北区江东南路24号
Address: No. 24 Tangmiao Rd., Hangzhou, Zhejiang

电话: 0571-85027738
Telephone: 0571-85027738
E-Mail: Locktest@126.com

浙江检验检疫局
@IDfKJ N . http://www.zjtj.cn:807
Online Business-Reception/Report Inquiries: http:// www.zjtj.cn:807

电话: 0571-85023552
Complaint Tel: 0571-85023552



Certificate of Compliance

No. 4Q200407M.WJPUT62

Certificate's
Holder:

Wenzhou Junyue Packing Co., Ltd.

The Building 5 in Yellow River Industrial Park, The
Century Avenue, Longgang Town, Cangnan County,
Zhejiang Province

Certification ECM
Mark:



Product:

Disposable Protective Mask

Model(s):

17.5 cm x 9.5 cm, 14.5 cm x 9 cm,
12 cm x 7 cm, KN95

Verification to:

Standard:

EN149:2001+A1:2009

related to CE Directive(s):

R 2016/425 (Personal Protective Equipment)

Remark: This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the technical documentation received from the manufacturer is satisfactory for the requirements of the ECM Certification Mark. The conformity mark above can be affixed on the products accordingly to the ECM regulation about its release and its use.

Additional information and clarification about the Marking:



The manufacturer is responsible for the CE Marking process, and if necessary, must refer to a Notified Body. This document has been issued on the basis of the regulation on ECM Voluntary Mark for the certification of products. RG01_ECM rev.3 available at: www.entecerma.it

Issuance date: 07 April 2020

Expiry date: 06 April 2025

Reviewer Technical
expert Amanda
Payne



Approver
ECM Service Director

Luca Bedonni





OFICINAS ESPAÑA

C/ Serrano 33, 4A28001 Madrid
C/ Poeta Querol 7, Pt. 13, 46002
Valencia

OFICINAS SUIZA

Theaterstrasse 5, 6003 Lucerne, LU

HORARIO EMAIL

Soporte 24x7 por el estado actual
causado por COVID-19
info@whitestoneeurope.com

www.whitestoneeurope.com